



STILLE imagiQ2™

Upgrades & accessories



Surgical perfection. For life.

imagiQ2™ standard configuration



Standard configuration

The STILLE imagiQ2™ comes with the following configuration:

- 300 kg / 661 lb weight capacity
- Ergonomic transport handles
- Battery life of up to 5 days*
- 4 swiveling wheels
- Pan handle
- Hand control + Emergency hand control
- Hand control extension cable
- Low attenuation mattress

**battery life contingent upon pt. weight and table use*

Part number	Description
530-1700	imagiQ2™ EU
530-1701	imagiQ2™ UK
530-1702	imagiQ2™ JP
530-1703	imagiQ2™ USA

imagiQ2™ upgrades to fit your needs

Mobility enhancement pack

When the table is used for frequent intra-clinic transportation, use the mobility enhancement pack for improved mobility. All wheels are designed to swivel and an additional fifth wheel is included to provide optimal mobility.

Configuration:

- 4 swiveling wheels and our proven linear 5th wheel.

Upgrade Ref: 535-1711

Load upgrade pack

Maximize the capabilities of the imagiQ2™ with increased patient weight capacity up to 661 lb / 300 kg. This premium upgrade has silver ergonomical transport handles.

Please note that this option must be specified on original order before delivery. It cannot be retrofitted to an already delivered unit.

Now included in standard configuration

Upgrade Ref: 535-1712

Performance pack

For maximum uptime, choose the performance pack. It includes a high-performance long-term battery pack with up to 5 days of durability in normal use.

This upgrade can be chosen at any time.

Now included in standard configuration

Upgrade Ref: 535-1710

Increased Trendelenburg pack

Maximize Trendelenburg capabilities by choosing the increased Trendelenburg pack. This allows you to extend the Trendelenburg angle from 20 degrees to 25 degrees. This upgrade is only compatible without the load upgrade pack.

Please note that this option must be specified on original order before delivery. It cannot be retrofitted to an already delivered unit.

Upgrade Ref: 535-1713

Upgrade features	Reference	Description
Mobility enhancement pack	535-1711	Extra 5th guiding wheel for intra-clinic transportation.
Increased Trendelenburg pack	535-1713	Increased Trendelenburg from 20 degrees to 25 degrees

imagiQ2™ accessory highlights

The imagiQ2™ offers a wide range of accessories for vascular, ERCP, pain management and electrophysiology; and it can be upgraded with features to adapt to surgical needs.

Products may not be available in all regions. Please contact your STILLE-representative for more information on market availability.



535-1765
Integrated control module
Combines the essential table control functions from the hand control and the pan handle into one device.

US-version: 535-1766.



535-1745
3D Head rest / extension
Carbon fiber head rest / extension that extends the radiolucent length with 30 cm (11.8”).



514-50-3SR040
Hand control holder
Holds hand control in 20° angle to provide one hand operation of the function buttons.

US-version: 514-50-3SR041.



535-1742
Adjustable procedure headrest
Fully radiolucent headrest, ensures easy intubation and superior carotid visualization. It fits both ends of the table and does not require side rails.



535-1721
Detachable side rail (80.5 cm)
Side rail enabling full table length side rails for accessory mounting. Not compatible with Adjustable head rest.

For US-version: 535-1720.

Also available in 70 cm (EU: 535-1729, US: 535-1723) and 28.5 cm (EU: 535-1725, US: 535-1724).



535-1750
Catheter tray 1,000 mm
Provides an extended work space, which is attached directly to the table.

Also available in 500 mm (535-1752) and 1,500 mm (535-1751).



535-1770
X-ray protective kit (grey)
Tailored and easy to fit protection against scattered radiation with lead rubber (0.5 Pb) material. For detailed information please see our brochure for Radiation Shields.

Also available in blue (535-1771).



535-1773
X-ray protective loop (blue)
Tailored and easy to fit protection against scattered radiation with lead rubber (0.5 Pb) material. For detailed information please see our brochure for Radiation Shields.

Also available in grey (535-1772) and as single panels to attach directly to the table top using Velcro strips (blue: 535-1774, grey: 535-1773)



535-1730
Radiolucent fistula arm board
Provides an unobstructed image and is connected directly to the carbon fiber top for Vascular procedures.



514-50-3SR025
Adjustable arm board
Nonradiolucent board provides ideal support for the patient's arm during IV.
US-version: 514-50-3SR027.



514-100530
Side gates, foldable (pair)
Side rail mounted gates to ensure the patient's safety and comfort during ERCP procedures. Recommended with side rails 535-1721 or 535-1720.
US-version: 514-100530-US.



514-100422
Widening sections (pair)
Side rail mounted side extensions, length 455 mm, with mattresses for lateral positioning of the patient. Other dimensions available upon request. Recommended with side rails 535-1721 or 535-1720.




514-100311
Anesthesia frame
 Enables hanging cloths or equal material to differentiate the sterile zone from the non-sterile zone.

US-version: 514-100311-US




514-100401
IV-pole
 Infusion pole with clamp for side rail fixation.

US-version: 514-100401-US



514-100554
Arm & leg holder
 Provides support to prevent the arm or leg of the patient from falling outside the operating table.



514-100552
Arm & leg holder curved
 Provides support to prevent the arm or leg of the patient from falling outside the operating table.



514-50-3SR032
Leg holders (pair)
 Leg holders Goepel type with clamps.

US-version: 514-50-3SR033.



505-1110
Dual arm board
 Radiolucent dual arm board for supporting the patient's both arms.



514-100465
Body/leg strap
 Strap to fix the patient, to be fixed to the side rails by Velcro.



514-3011-00
DORO® Side Rail Adaptor Amsco
 Side rail mounted adaptor for DORO Headrest (non-radiolucent).

US-version: 514-3011-10



514-100334
Head cushion
 Cushion to stabilize the patient's head during procedures.



535-1762
Foot control
 Compact multipedal footswitch for discrete placement beneath the table with soft touch buttons and a light operating force.



514-50-3SR036
Adjustable tray with clamp
 Adjustable tray to place instruments and other important devices useful to the surgical team

US-version: 514-50-3SR37



535-1760
Extra pan handle
 Extra pan handle allows panning the table from different sides.

US-version: 535-1761

Additional accessories

Item	Part number (EU)	Part number (US)
Head rest, non-adjustable	535-1741	535-1741
Lateral support with clamp	514-50-3SR028	514-50-3SR030
General purpose equipment drapes (50 pcs)	514-17500	x
Basic arm board	x	512-5000



STILLE imagiQ2™ | Accessories

STILLE – Surgical perfection

STILLE is one of the oldest medical device companies in the world. At STILLE, we develop, produce and sell high-quality surgical instruments and surgical imaging tables to physicians worldwide, and has done so for 180 years.

STILLE's offering consists of a strong brand and products renowned for their superior quality. We manage our business in two distinct business areas – Surgical Instruments and Surgical imaging tables.

Surgical instruments

In our Surgical Instruments business area, we have a complete offering of handcrafted, forged instrumentation for all fields of open surgery, with a special emphasis on cardiovascular surgery, and reconstructive and aesthetic plastic surgery.

Our instruments feature delicate design, industry-leading durability and a unique feel, which has made STILLE a well-known and highly respected supplier by surgeons throughout the world.

Surgical imaging tables

In our Surgical imaging tables business area, we offer specialty tables for minimally invasive procedures, with pricing, safety and features optimized for both in- and outpatient settings. The STILLE imagiQ2™ surgical imaging table has the highest radiolucency in the world, and is optimized for advanced mobile Hybrid OR settings as well as peripheral and interventional procedures.

The STILLE Medstone line of imaging tables offer reliability and optimal performance at an affordable cost with the highest standards of Swedish technology.

Learn more at www.stille.se



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30
Torshälla, Sweden
Phone: +46 8 588 580 00
Fax: +46 8 588 580 05
info@stille.se

www.stille.se

STILLE Surgical Inc.

840 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.665.1614

Distributed by: